

Estimation of the total amount of Ca^{2+}

At any pH beyond 10, Mg^{2+} gets precipitated as $\text{Mg}(\text{OH})_2$. So, the solution contains only Ca^{2+} after making the pH over 10 by the addition of NaOH . Then the Ca^{2+} present in the resultant solution can be estimated by complexometric titration with EDTA solution.

Requirements:

Water sample, EDTA solution (0.01M), Eriochrome Black T, Buffer solution, NaOH solution, Burette, Pipette, conical flask (100 mL), Beaker (200 mL).

Procedure:

a) Standard 0.01 M EDTA Solution was prepared by weighing about 3.8 g of the disodium EDTA salt ($\text{Na}_2\text{H}_2\text{Y} \cdot 2\text{H}_2\text{O}$) into a 1-liter volumetric flask followed by its dissolution and dilution to the mark with deionized water.

c) Estimation of the total amount of Ca^{2+} & Mg^{2+}

The burette was filled with standard EDTA solution to the zero level, following usual precautions. 10 mL of the given water sample is pipetted out into a clean conical flask. 2-3 mL ammonia buffer and 2 drops of EBT indicator are added and titrated against EDTA from the burette. The endpoint was the change of color from wine red to steel blue. The titration is repeated to get three concordant titer values.

d) Estimation of the total amount of Ca^{2+}

In a conical flask or beaker 10 mL water sample was taken and pH was made beyond 10 by the dropwise addition of NaOH solution. The solution was shaken and filtered. The filtrate was titrated by EDTA solution using EBT as the indicator.

Results and calculation:

Titration-1 Estimation of Ca^{2+} & Mg^{2+} by titration with EDTA

The volume of the given water sample (mL)	Burette Reading (mL)			The volume of EDTA solution (mL)
	Initial	Final	Use volume	
10	0	0	12	12
10	12	24.1	12.1	12.1
10	24.1	36.1	12	12

$$\text{Mean} = \frac{12 + 12.1 + 12}{3} = 12.033$$

Titration-2 Estimation of Ca^{2+} by titration with EDTA

The volume of the given water sample (mL)	Burette Reading (mL)			The volume of EDTA solution (mL)
	Initial	Final	Use volume	
10	0	9	9	9
10	9	18	9	9
10	18	27.1	9.1	9.1

$$\text{Mean} = \frac{9 + 9 + 9.1}{3} = 9.033$$

Calculation:

1 ml of 0.01 M EDTA \equiv 1 mg of CaCO_3

9.03 ml of EDTA \equiv 9.03 mg of CaCO_3

Calculation of the amount of Ca^{2+} & Mg^{2+}

Volume of sample water taken = 10 ml

Volume of EDTA solution consumed = 12.033 ml

$$\begin{aligned} \text{Total amount of } \text{Ca}^{2+} \text{ \& } \text{Mg}^{2+} &= \frac{\text{Volume of EDTA solution consumed} \times 1000}{\text{The volume of the hard water taken}} \text{ ppm} \\ &= \frac{12.033 \times 1000}{10} \text{ ppm} = 1203.3 \text{ ppm} \end{aligned}$$

$$\begin{aligned} \text{Total amount of } \text{Ca}^{2+} &= \frac{\text{Volume of EDTA solution consumed} \times 1000}{\text{The volume of the hard water taken}} \text{ ppm} \\ &= \frac{9.033 \times 1000}{10} \text{ ppm} = 903.3 \text{ ppm} \end{aligned}$$

$$\begin{aligned} \text{Total amount of Mg}^{2+} &= (\text{Total amount of Ca}^{2+} \text{ \& Mg}^{2+} - \text{Total amount of Ca}^{2+}) \\ &= \underline{300} \text{ ppm} \end{aligned}$$

Result:


The collected water sample contains

$$\text{Total amount of Ca}^{2+} \text{ \& Mg}^{2+} = 1203.3 \text{ ppm}$$

$$\text{The total amount of Ca}^{2+} = 903.3 \text{ ppm}$$

$$\text{The total amount of Mg}^{2+} = 300 \text{ ppm}$$

[Signature]
19/06/23



EBT

~~NOT~~

EDTA

